

## **Aristo M1457-3 Series**

Double Ring Microscope Lamp  
for Bausch & Lomb and Olympus Microscope

Specifications:

- Reflector: 4"
- ID: Bausch & Lomb: 1.5" or  
Olympus: 1.89"
- Connection to power pack: 3-pin
- Power pack to outlet:  
3-pin grounded connection



Complete System shown:  
Reflector, Lamp, Power Supply, Power Cord with 3-pin black connector  
and either 1.5" or 1.89" adapter



Replacement Lamp  
Head Assemble shown:  
Reflector, Lamp, Power  
Cord and Connector



Replacement Lamp:  
Available in W31, W45  
and W55 colors



Replacement Power Pack shown:

### **Caution:**

*Instructions are for guidance only. Aristo/Voltarc will not be responsible for any damage to lamp or lamp housing for items improperly installed.*





1. Cut electrode lead wires to 1/2 " long and fold back against the electrode.
2. Unscrew electrode shield from lamp housing and pry the lamp out of it's fuse clips.
3. Slide the pressure fittings over the electrodes making good contact.
4. Carefully position the lamp over the fuse clips and ease them back into the fuse clip.

